



Be Right™



EDTA Digital Titrator Cartridge, 0.800 M

Product #: 1439901
CAD Price: \$42.30
Available

For determination of Calcium and Total Hardness by Digital Titrator titration.

0.800 M EDTA. Digits = mg/L CaCO₃ in 100 mL sample.

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	EDTA (ManVer 2)
Method Name:	EDTA (CalVer 2)
Method Number 2:	8213
Method Number 3:	8329
Number of tests:	40 - 100
Parameter:	Durezza, total, secuencial
Platform :	Digital Titrator
Range:	100 - 4000 mg/L CaCO ₃
Titrant:	0.800 M EDTA

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)