



Boson COVID-19 Test Kit - 20 Pack / White



Using a nasal swab to get sample fluid, antigen tests can produce results in minutes. Because these tests are faster and less expensive than molecular tests are, experts consider antigen tests more practical to use for large numbers of people. Since the antigen test can be high specificity, a positive antigen test result is considered very accurate. Kit Contents: 20 Rapid SARS-CoV-2 Antigen Test Cards (Individually packaged), 20 Buffer Packs, 20 Sterilized nasopharyngeal swabs (Individually packaged), 20 Extraction tubes (Individually packaged), 1 Instruction for use, 1 Tube rack (for 6 tubes/rack)

Key Features

- Sample type: Nasopharyngeal swab
- Storage temperature: 2–30°C
- Shelf life: 18 months
- Test time: 15–20 minutes (do not read results after 20 minutes)

| | | | |
|----------------------|---------------------|---------------------|---------------------|
| SKU | 1N40C5 | Inner Dims IN | 8"L x 4.9"W x 2.5"H |
| Master Dims IN | 8"L x 4.9"W x 2.5"H | Inner Weight Lbs IN | 0.4 |
| Master Weight Lbs IN | 0.4 | Inner Pack QTY | 20 Kits |
| Master Cubed IN | 98.0 | Inner UPC | 837528001437.0 |
| Master Cubed FT | 0.06 | Master Case Units | 20 Kits |
| Master Case QTY | 1 Box | | |
| Master UPC | 10837528001434 | | |
| Master Skid Qty | 1250 | | |
| Master Tie | 250 | | |
| Master High | 5 | | |

Available Docs: MDEL LIC, MDEL LIC II

Right Product. Right Price.  